


- * Tubular mattress is more suitable & effective for patients having bed sores from stage 1-3
- * It has a better support & can handle patients having weight upto 125 kg.
- * It is orthopaedically designed for better body comfort.
- * Deeper air cells with better support for body pressure. More durable & longer life.
- * Individual cells can be replaced instead of complete mattress.



Specifications

Air Pump

Power :220 v/50 Hz

Pressure Range :30 mm Hg -70 mm Hg

Air Output :4 -4.5 litre/minute

Cycle Time :10 minutes

Suspension system design for motor .

No noise and vibration.

PVC Mattress

- Made of heavy gauge, high elastic, low toxic Nylon/PVC
- Reusable and tear-resitant
- Dimensions : 190cm Lx90cm Wx 10cm (inflated)
- Dimension of individual Tube = 90 x 10cm
- No. of tubes in one mattress = 18