



Ready baby Manual Breast Pump

with Silicon Breast Shield

New



Brand-new design of Ergonomics pull handle makes the suction more powerful, and let you feel no any hand tingle

Mother's Milk is like Nectar

Give exclusive breastfeed to newborn up to six months

Remember, Only Mother's Milk gives the baby complete food and protects from diseases

Mother must give her first thick, yellow milk to the newborn immediately after birth

Breastfeeding should start within one hour of birth





Ready baby Manual Breast Pump

For Hospital, Clinics and Home Use

Function & Features:

Breast feeding is the best choice for babies' growth. Newborn can obtain immunogenic molecules from breast milk, but sometimes you find it is inconvenient to feed your baby directly. Our breast pump is the ideal appliance.

The manual breast pump adopts pull handle design, users can make breast sucking at their comfortable gesture; its dual-ring silica gel and unique 4 air inlet design makes suction around mamilla uniformly and alleviate pain and avoid the mamilla collapses, let mothers healthier and more beautiful.



List Of Components

1. Lever Handle
 2. Connector
 3. Silicon Adaptor shield
 4. Funnel
 5. 2-Component Connector
 6. Value Head
 7. Value Membrane
 8. Feeding Bottle
 9. Bottle Cover
 10. Nipple
 11. Bottle Connector
 12. Bottle Stand
- Accept 1 can not be sterilized, the rest all can be cleaned and sterilized in the water.
- When using feeding bottle sterilizer and dryer, please read the specifications or putting them into hot water for sterilizing in 3 minutes.

Thank you for using our products, we are committed to offer you better quality, much more considerable service and favorable prices to consumers.



Exclusively Imported & Mkt. By :
Mehar Healthcare Corporation

Admn. Off: 16/89, Lane No. 4, Faiz Road, Karol Bagh, New Delhi - 110005
Phone : 011-23520753, 23539529, 23864115, 23861140
Fax : 011-25812657 Mobile : 98101-00623
E-mail : meharhealthcare@rediffmail.com
www.meharhealthcare.com

New



Gentle Reminder

- Put it in the place where the children can not reach it.
- Please sterilize all the components before use.
- Moms sit comfortably to make breast sucking conveniently.
- Please cover with the dust-proof cover after use.
- Suggestion: Please clean and sterilize the product after use every time.

Product Features

Manual Breast Pump is designed to help mothers express breast milk with ease. It is light and portable, allowing mothers to use it conveniently be it indoor or outdoor. Feeding bottle, bottle cover set and teat are included to make it a practical and complete breast pump set. A bottle holding base is added into the set to prevent the feeding bottle from tipping thus no milk spillage. Another safe, reliable and durable product only from baby and mother care product ranges.

Technique Parameter

Stimulation phase: 0—150mmHg
0—20KPa
Expression phase: 0—230mmHg
0—20KPa



small and exquisite
convenient for carrying

massage function

uniform suction

pull handle

