

Digital Blood Pressure Monitor

Wrist/ Full automatic



CH-657

Measurement system	Oscillometric method
Measurement range	Pressure : 0 ~ 280 mmHg Pulse : 40 ~ 180 pulse/min
Accuracy	Pressure : ± 3 mmHg Pulse : ± 5 %
Power source	AAA size $\times 2$
Battery life	Alkaline battery : Approx. 400 times
Rated voltage	DC3V \approx 3W (\approx : Direct Current)
Outside dimensions	W69 \times H70 \times D15 mm
Weight	Approx. 86g (without batteries)

- Irregular Heartbeat Indicator
- Body movement indicator
- 90 memories function
- Average of last 3 readings
- Easy to read LCD display
- Automatic power off function
- Handy and slim design
- With compact storage case
- Simultaneous blood pressure and pulse reading

