

New



CH-456

CE 0086

- Irregular Heartbeat indicator
- Body movement indicator
- 90 memories function
- Average of last 3 readings
- Easy to read LCD display
- Cuff storage holder
- Automatic power off function
- Simultaneous blood pressure and pulse reading
- Time and date
- AC adapter is optional (AC-230CZ)

Measurement system:	Occlusometric method
Measurement range:	Pressure: 0-280 mmHg Pulse: 40-180 pulse/min
Accuracy:	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Power source:	AA size x 4 Optional AC/DC adapter (DC 6V 0.6A)
Battery life:	Manganese battery / Approx. 150 times Alkaline battery / Approx. 300 times
Rated voltage:	DC5V (1.5V x 3) (Direct Current)
Outside dimensions:	W125 x H146 x D122 (mm)
Weight:	Unit / Approx. 380g (without batteries) Cuff / Approx. 130g



CH-432

CE 0086

- 90 memories function
- One touch fully automatic
- Average of last 3 readings
- Low Battery consumption
- Silent Inflation

SPECIFICATIONS

Measurement System:	Occlusometric method
Measurement range:	Pressure: 0-280 mmHg Pulse: 40-180 pulse/min
Accuracy:	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Power Source:	AA Size x 4
Battery life:	Manganese battery: Approx. 300 times Alkaline battery: Approx. 1000 times
Rated Voltage:	DC5V (1.5V x 3) (Direct Current)
Outside dimensions:	129mm(W) x 55mm(H) x 119mm(D)
Weight:	Unit: Approx. 245g (without batteries) Cuff: Approx. 130g



Your Area Distributor

CH-456 **New**

CITIZEN
Micro HumanTech

Digital Blood Pressure Monitors

Monitor Your Blood Pressure Regularly,
Hyper Tension is a Silent Killer



CE 0086

How to operate



1 Wrap the Cuff around your arm and make sure the 'blue' marker is located in the center of your arm.

2 Place your arm on a table with your palm opened and relaxed. Press the START/STOP switch to START measurement.

3 The blood pressure generally varies in a day. Preferably, you should measure your blood pressure at the same time each day wherever possible.



Exclusively Imported & Serviced By:

Mehar Healthcare Corporation

Admin. Off: 18/50, Lane No. 4, F-4 Road, Karol Bagh, New Delhi - 110028

Phone: 011-25820753, 25538929, 23864115, 23861140

Fax: 011-25812657, Mobile: 98101-06623

E-Mail: meharhealthcare@rediffmail.com, meharhealthcare@gmail.com

Website: www.meharhealthcare.com