

Ear type

new



CE
0197

CTD504 CELSIUS (°C) · FAHRENHEIT (°F)

Measurement range	34.0°C ~ 43.0°C (93.2°F ~ 109.4°F)
Accuracy	Celsius: ± 0.2°C (35.5 ~ 42.0°C) Fahrenheit: ± 0.4°F (95.9 ~ 107.6°F) Other range: ± 0.3°C (± 0.5°F)
Display	3digit with decimal point . °C, °F, Low Battery Symbol mark, Too High Ambient Temperature mark, Too Low Ambient Temperature mark, Ready to Use Symbol mark, Memory Mode Symbol mark.
Battery	DC3V (Lithium Battery : CR2032 × 1)
Battery life	1500times
Outside dimensions	Approx. W39 × L120 × 52mm
Weight	Approx. 68.6g (Including battery)

- Antibacterial Ear Thermometer
- No need to use probe cover
- 1 ~ 3 second measurement
- 12 memory
- °C / °F switching
- Auto power off




CAUTION

Measurement result with this unit may be different from a result by other type of thermometers due to a different measuring place.