

Function mark explanation

Celsius/ Fahrenheit



Celsius



Fahrenheit



°C / °F
Conversion

Measurement Method Measurement Point



Earlobe
measuring

Main Features and Specifications



Display



Water
resistant



Anti bacterial
material



Auto
power off

Memory and Notification



12 memory



Memory
recall



Buzzer
notice



Probe cover is
not necessary



Measuring
Speed

Battery type Battery Life (Approximate Value)



Battery life



Accessories



Case



Stick type

new



CTA301 CELSIUS (°C) · FAHRENHEIT (°F)

Measurement method Temperature peak hold thermometer
Measurement range 32.0°C ~ 42.9°C (90.0°F ~ 109.9°F)
Accuracy Celsius : ± 0.1°C (35.5 ~ 42.0°C)
Fahrenheit : ± 0.2°F (95.9 ~ 107.6°F)
Display 3 digits with decimal point, °C, °F,
battery weak mark and memory mark
Battery DC1.5V (Micro Alkaline Battery : LR41 × 1)
Battery life Approx. 2 years (If used 10 minutes a day)
Outside dimensions Approx. W19 × L125.7 × D11.2 mm
Weight Approx. 9.7g (including battery)



- °C / °F switching
- General purpose use
- Auto memory recall of last reading
- Buzzer notice
- 3 digits display
- Auto power off



new



CTA302 CELSIUS (°C) · FAHRENHEIT (°F)

Measurement method Temperature peak hold thermometer
Measurement range 32.0°C ~ 42.9°C (90.0°F ~ 109.9°F)
Accuracy Celsius : ± 0.1°C (35.5 ~ 42.0°C)
Fahrenheit : ± 0.2°F (95.9 ~ 107.6°F)
Display 3 digits with decimal point, °C, °F,
battery weak mark and memory mark
Battery DC1.5V (Micro Alkaline Battery : LR41 × 1)
Battery life Approx. 2 years (If used 10 minutes a day)
Outside dimensions Approx. W19 × L126.7 × D11.2 mm
Weight Approx. 9.7g (including battery)



- °C / °F switching
- General purpose use
- Water resistant & fully washable
- Auto memory recall of last reading
- Buzzer notice
- 3 digits display
- Auto power off



new



CTA303 CELSIUS (°C) · FAHRENHEIT (°F)

Measurement method Temperature peak hold thermometer
Measurement range 32.0°C ~ 42.9°C (90.0°F ~ 109.9°F)
Accuracy Celsius : ± 0.1°C (35.5 ~ 42.0°C)
Fahrenheit : ± 0.2°F (95.9 ~ 107.6°F)
Display 3 digits with decimal point, °C, °F,
battery weak mark and memory mark
Battery DC1.5V (Micro Alkaline Battery : LR41 × 1)
Battery life Approx. 2 years (If used 10 minutes a day)
Outside dimensions Approx. W19 × L126.7 × D11.2 mm
Weight Approx. 10.2g (including battery)



- °C / °F switching
- General purpose use
- Water resistant & fully washable
- Auto memory recall of last reading
- Buzzer notice
- 3 digits display
- Auto power off
- Antibacterial Thermometer

