

# Nebulizer

## Ultrasonic type



### CUN60-CE CUN60-120V

|                         |   |
|-------------------------|---|
| Power                   | AC220-240V, 50/60Hz<br>AC100V, 50/60Hz  |
| Power consumption       | Approx. 10W   |
| Ultrasonic Frequency    | 2.5 MHz   |
| Nebulization rate       | 3-speed adjustable  |
| Particle size           | 5 $\mu$ m MMAD  |
| Medication cup capacity | Max. 8ml  |
| Auto shut-off timer     | Approx. 10min   |
| Dimensions              | W88 x H143 x D53 mm   |
| Weight                  | Approx. 350g<br>( Including rechargeable battery-pack )   |
| Accessories             | Nose piece, Mouth nozzle, Extension tube,<br>Adult & child masks, AC adapter,<br>Medication cup ( 5 pcs. ), Soft carrying case. |

- Handy design  
with built-in detachable rechargeable battery pack,  
compact and light-weight
- Easy to use  
uses water as ultrasonic conducting medium, no gel is  
required
- Silent operation  
thanks to the ultrasonic vibration technology
- 3-speed nebulization rate control  
provides optimal treatment for different density
- Built-in timer  
automatically shuts off after 10 minutes use



Totally Portable