



Ready Thermo™ Flexy

Digital Thermometer

FLEXIBLE TIP
WATER PROOF



DUAL
SCALE C/F

JAPANESE
TEMPERATURE
SENSOR

BEEP AT PEAK
TEMPERATURE

EASY TO READ
LCD DISPLAY

A high quality healthcare product
for your whole family





Ready Thermo™ Flexy

Digital Thermometer

- WATER PROOF • FLEXIBLE TIP
- DUAL SCALE C/F • BEEP AT PEAK TEMPERATURE
- ON / OFF BUTTON AUTO SHUT-OFF
- LAST TEMPERATURE MEMORY RECALL
- FAST ACCURATE READOUT IN 30 SECOND
- EASY TO READ LCD DISPLAY
- CONVENIENT CLEAR STORAGE CASE
- JAPANESE TEMPERATURE SENSOR
- MEETS ASTM STANDARD E 1112
- THE TEMPERATURE SENSOR USED IN THIS THERMOMETER HAS BEEN MADE IN JAPAN



Specifications

Measurement range	90.0°F (32.0°C) ~ 109.9°F (42.9°C)
Display resolution	0.1°F (0.1°C)
Accuracy	+/-0.2°F.98.0°F-102.0°F (+/-0.1°C, 37.0°C-39.0°C) +/-0.3°F.96.4°F-98.0°F(+/-0.2°C, 35.8°C-37.0°C) +/-0.3°F. 102.0°F-106.0°F (+/-0.2°C, 39.0°C-41.0°C)
Average measurement time	Oral: 30 seconds; Rectal: 30 seconds; Axillary: 60 seconds.
Auto power off	10 minutes ± 2 minutes
Battery	LR41 or SR41
Dimension	133*21*11mm
Operating condition	Temperature: 50.0°F ~ 104.0°F (10°C ~ 40°C) Humidity: 30% ~ 85% RH
Storage condition	Temperature: 14.0°F to 140.0°F (-10°C ~ 60°C) Humidity: 25% ~ 90% RH

Exclusively Imported By :



Mehar Healthcare Corporation

Admn. Off :16/89, Lane No. 4, Faiz Road,
Karol Bagh, New Delhi - 110005
Phone : 011-23520753, 23539529, 23864115, 23861140
Fax : 011-25812657 Mobile : 98101-00623
E-mail : meharhealthcare@rediffmail.com
www.meharhealthcare.com

Distributor

