



Anti Decubitus Alternating Air pressure pump
& mattress system (Bed Sore Cure System)

Ready Bed

ECO



- Quiet Pump
- Low Vibration
- Low Temperature



Anti Decubitus Alternating Air pressure pump
& mattress system (Bed Sore Cure System)

Ready Bed

ECO

mattress:

bubbles on A-B lines alternating
inflating provides
superior therapeutic treatment
made of durable heavy,
low toxic premium vinyl

Mattress Characteristics

Max. User Weight : 90-100 Kgs
Mattress Dimensions : 200*90cm
Materials : PVC 0.3mm
Reusable and tear-resistant
With 130 bubble cells.



Pump:

Simplicity, reliability and
easy of use
Aluminium motor cover
quiet pump low noise,
low vibration, low temperature
with knob easy adjust pressure
fixed holk easy fix on
bed frame

Compressor Characteristics

Cycle Time : 6 Mins.
Air Output : 6 - 7 L / Min.
Pressure Range : 35 - 70 mmHg
Electrical Power Supply : 220 v / 50 Hz
Protection Against Class II
Electric Shocks



Exclusively Imported & Mkt. By:
Mehar Healthcare Corporation

Admn. Off: 16/89, Lane No. 4, Faiz Road, Karol Bagh, New Delhi - 110005
Phone: 011-23520753, 23539629, 23864115, 23861140
Fax: 011-25412957 Mobile: 98101-00623
E-mail: meharhealthcare@rediffmail.com
www.meharhealthcare.com

DEALER Stamp :

