

ULTRASONIC NEBULIZER



Ready Flow Heavy Duty Ultrasonic Nebulizer (For Clinic & Hospital)

- Totally Noiseless
- Continuous Operation \geq 4 Hours
- Large Medicine Capacity

Main Technical Parameters

- | | |
|---|--|
| • Ultrasonic frequency | 1.7 MHz \pm 10% |
| • Continuous working time | ≥ 4 hours |
| • Consumed power | 4.90 VA |
| • Maximum capacity of small atomizing cup | 150 mL |
| • Maximum capacity of big atomizing cup | 300 mL |
| • Maximum atomizing rate | 3.4mL / min |
| • Power supply voltage | ~ 220V, 50 Hz |
| • Toning range | stepless adjustable 0 - 80 min |
| • Water tank capacity | 300 mL (submerged by water level line) |
| • Particle Size | 0.5 to 5 μ m |



New Totally portable!

Ultrasonic Nebulizer

Handy design and light in weight with rechargeable battery pack

- Handy design... with built-in detachable rechargeable battery pack, compact and light-weight, users water as ultrasonic conducting medium, no gel is required thanks to the ultrasonic vibration technology.
- Easy to use...
- Silent operation...
- 3-speed nebulization rate control... provides optimal treatment for all kinds of medication with different density, automatically shuts off after 10 minutes
- Built-in timer...
- Electromagnetic compatibility... complies with EN60601-1-2 requirement as specified in EC directive 93/35/EEC

SPECIFICATIONS

Power
Power consumption
Ultrasonic Frequency
Nebulization rate
Particle size
Capacity
Auto shut-off timer
Dimension
Weight
Accessories

AC220-240V, 50/60Hz
Approx. 10W
1.7MHz
3-speed adjustable
SUNMAD Medication cup
50x 50x
Approx. 1.5kg
300x200x100mm (DxLxH)
Approx. 300g (including rechargeable battery pack)
Nose piece, Mouth-piece, Exhalation tube, Adult & Child masks, AC adapter, Medication cup (Spray), Soft carrying case.



0434



Mehar Healthcare Corporation

Your Most Reliable Source for
NEBULIZERS

Introducing

Ready Flow

Heavy Duty Ultrasonic Nebulizer
(For Clinic & Hospital)

Totally Noiseless

Continuous Operation \geq 4 Hours

Large Medicine Capacity



ReadyJet™
NEBULIZER CHAMBER



Imported & Marketed By :

Mehar Healthcare Corporation

Admn. Off-16/89, Lane No. 4, Faiz Road Karol Bagh, N Delhi-5

Phone : 011-23520753, 23539529, 23864115, 23861140

Fax : 011-25812657 Mobile : 98101-00623

E-mail : meharhealthcare@rediffmail.com • www.meharhealthcare.com