

Mehar Healthcare Corporation



# Ready Hospy

High Quality, Heavy Duty, Powerful  
Piston Compressor Nebulizer  
for Hospital, Clinics & Professional Use



- With Tough Metal Body
- Adjustable Output Pressure
- With Pressure Display Meter
- For Long Term Continuous Use
- With Thermal Protector

#### TECHNICALS FEATURES :

▶ Power	:	220V 50 Hz 180 VA
▶ Particle Size	:	0.5 to 10 $\mu$ m
▶ Nebulization Rate	:	$\geq$ 0.2ml/min
▶ Liter Flow Range	:	5.0 to 8.0 Lpm
▶ Operating Pressure Range	:	14.5 to 25.0 PSI (0.1 to 0.17Mpa)

CE  
0434